Imperial College London


Design and Professional Practice 2

Risk Analysis


Dr Ian Radcliffe
























Imperial College London _{S-Sever}


S - Severity


4	Rating S	Criteria: Severity of effect	Consequence	Treatment			
	10	Death	-	-			
	9	Quadriplegia	Life-long medical care necessary / coma / permanent damage				
	8	Amputations, paraplegia, blindness, deafness, traumatic brain injury (severe), fourth-degree burns	Life-long medical care necessary / coma / permanent damage	Hospital stay			
	7	Complex fractures, open fracture, inner injuries, traumatic brain injury (severe), third-degree burns	Permanent damage possible	Hospital stay			
	6	Gash, fractures, torn muscles, articular cartilage injury, traumatic brain injury (moderate), second-degree burns	Permanent damage possible	Hospital stay			
	5	Gash, fractures, torn muscles, articular cartilage injury, traumatic brain injury (mild), second-degree burns	Reversible injury	Hospital stay or ambulant treatment			
	4	Severe cuts, severe scratches, severe contusions, strains, first-degree burns	Reversible injury	Ambulant treatment or self-treatment			
	3	Minor cuts, minor scratches, minor contusions, stiff muscles, tension, blisters, excoriations, sickness, first-degree burns	Discomfort during application up to three days after application	Self-treatment			
	2	Slight sickness, pressure marks	Discomfort	-			
	1	No harm	-	-			

Risk Evaluation









Imperial College London

Risk Management

	Risks	s Effects / Consequences S1 O1 D1 RPN Preventing measures							00	Da	
8	RISKS	Effects / Consequences	S1	OI	וט	before	Preventing measures	S2	02	DZ	RPN after
	Operational hazards										
	Wrong placement of connections between the device and the oxygen supply	Lack of gas delivered to the baby and hence possible asphyxia and death.	10	3	3	90	Introduce labels at the connections to ensure complete clarity and distinction between input and output.	10	2	2	40
	Incorrect placement of oxygen mask or nasal cannulae from the controller to the neonate	Lack of gas delivered to the baby and hence possible asphyxia and death. Also, possible physical damage of the neonate.	10	3	3	90	Add an extra section in the user manual to explain how to place the nasal cannulae or oxygen mask to the baby.	10	2	2	40
	Ambient light interference	Incorrect measurements from pulse oximeter.	10	6	8	480	Include machine learning in the code used with the pulse oximeter to detect unusual behaviour of the sensor. Include a warning in the user manual advising staff to cover the pulse oximeter from ambient light.	10	6	1	60
	Movement of pulse oximeter	Generates incorrect measurements and could potentially block the valve's code.	10	9	4	360	Improve the design of the strap connecting the pulse oximeter to the neonate. Include machine learning in the code used with the pulse oximeter to detect unusual behavior of the sensor.	10	6	1	60

